

MICROMETER Pneumatic & Manual

Code. N1102 and N1102P

Specially suitable to measure the thickness of various materials, paper, tissue, etc.

Aplicable Standard : ISO 534, ISO 12625-3, TAPPI T411



Foot switch for pneumatic model

PTA SPAIN

C/ Gurutzegi N° 16 Pol. Belartza
20018 SAN SEBASTIÁN
(Guipúzcoa) - Spain

PTA FRANCE

14, rue Evariste Galois / ZAC les Ruieres
38320 EYBENS - France

PTA ASIA

50 East Coast Road #02-73 Roxy Square
428769 SINGAPORE - Singapore

**DEVICE DESCRIPTION :**

The robust micrometer is available in either pneumatic or manual configuration. Depending on the model, the device is operated via a hand lever or a pneumatic foot pedal. The descending speed acc. to standard is set with a throttle. A battery powered digital measurement gauge is mounted on the device for displaying the measurement results.

TEST DESCRIPTION :

Depending on configuration, the measuring pin is lifted via the pneumatic foot pedal or manually with the lever mounted on the side, and the sample is placed in the measurement area. Then the pedal or lever is released, and the measuring pin lowers at a speed limited by the throttle onto the sample. The digital measurement gauge displays the result of the measurement.

SPECIFICATIONS :

- ✓ Bound the lowering speed via a throttle
- ✓ Battery operated digital measurement gauge
- ✓ Resolution: 1 μm , reproducibility : $\pm 2 \mu\text{m}$

MODELS :

N 1102 - Manual model : 10 cm² measuring pin

N 1102P - Pneumatic model : 10 cm² measuring pin

PTA SPAIN

C/ Gurutzegi N° 16 Pol. Belartza
20018 SAN SEBASTIÁN
(Guipúzcoa) - Spain

PTA FRANCE

14, rue Evariste Galois / ZAC les Ruieres
38320 EYBENS - France

PTA ASIA

50 East Coast Road #02-73 Roxy Square
428769 SINGAPORE - Singapore



PAPER TESTING ASSOCIATION

CORRUGATED BOARD Testing Equipment



CONNECTIONS :

Electricity : battery

Compressed air : approx. 400 kPa (only on the pneumatic model)

WEIGHT AND DIMENSIONS (ONLY EQUIPMENT) :

200 x 250 x 300 mm (WxDxH) / 12 Kg

WEIGHT AND DIMENSIONS WITH PACKING BOX :

400 x 400 x 400 mm (WxDxH) / 16 Kg

PTA SPAIN

C/ Gurutzegi N° 16 Pol. Belartza
20018 SAN SEBASTIÁN
(Guipúzcoa) - Spain

PTA FRANCE

14, rue Evariste Galois / ZAC les Ruieres
38320 EYBENS - France

PTA ASIA

50 East Coast Road #02-73 Roxy Square
428769 SINGAPORE - Singapore

info@ptagroup.com

www.ptagroup.com