




NOVILOG | Log-03

Mesure n°	Force (N)	Temps (s)	Temps (s)
Mesure n° 1	4.801	14.875	14.875
Mesure n° 2	4.801	14.875	14.875
Mesure n° 3	4.801	14.875	14.875
Mesure n° 4	4.801	14.875	14.875
Mesure n° 5	4.801	14.875	14.875
Mesure n° 6	4.801	14.875	14.875
Mesure n° 7	4.801	14.875	14.875
Mesure n° 8	4.801	14.875	14.875
Mesure n° 9	4.801	14.875	14.875
Mesure n° 10	4.801	14.875	14.875

Description

Novilog is a software measurement for your all devices.

Novilog is a quality control software for laboratories. Its operation allows it to be adapted to many situations:

- it can be used as a simple acquisition software for one device,
- it can be linked to a large number of devices on a single PC,
- it can run in cloud mode and connect over the network with as many stations and devices as desired.

NOVILOG | Log-03

Technical specifications

Main properties

- Compliance of measurements according to ISO.
- Security and traceability of measurements according to ISO.
- Centralization of data according to ISO.
- Production control.
- Calculations and statistics.
- Connection to all devices in our range.
- Connection of other devices possible on condition that we can be provided with the connection protocol of the devices to be connected. Extra charge possible on quotation.
- Multi-language software.

Features

- Recovery of all the measurement data of your devices.
- Exploitation of results with adapted statistics and exports.
- Interconnection to other systems to exchange data via different protocols.

Characteristics

Minimum resolution 1280x 800

Contents of the delivery

A USB key
A dongle